



Indi™

The power to express myself



Indi is the world's first all-in-one, consumer-focused speech tablet for symbol-based communicators. Unlike anything else on the market, Indi was designed specifically for Augmentative and Alternative Communication (AAC) and comes with all of your needs fully integrated, like our intuitive and easy-to-use Snap™ software. At the heart of Snap is the Core First® pageset, Tobii Dynavox's systematic delivery of a Core Word framework that continuously grows vocabulary and literacy skills, while empowering communication success and independence.



Built for your world

with an easy-to-hold, ergonomic design, Indi is ultra-portable and ready to go out-of-the-box.



Adapts to your access method

with numerous compatible inputs such as switch-scanning, head-mouse, and touch with keyguards, etc.

Indi is designed around your needs



Created to be heard

with powerful, outward facing, built-in speakers providing clear speech, exceptionally loud output, and a wide range of authentic sounding voices.



Made for everyday interactions

with capabilities that go beyond speech communication like Environmental Control Units (ECUs), Infrared (IR), access to social media, email, and texting.

Tobii Dynavox gives you the power to focus on what matters

Our goal is to help you focus on communication growth, and not worry about technology. To do this, Indi users receive access to a truly unparalleled ecosystem of content and a library of resources. Free resources include Pathways for Core First, Discover Tobii Dynavox, Core First Learning, Core First Books & Lessons, Boardmaker® Student Center, and more.

Indi with Snap + Core First, all for \$999

Indi, our consumer-focused (non-funded) speech tablet includes everything needed to start right out of the box and build success, including our Snap + Core First communication software. Ordering Indi is as easy and fast as purchasing a consumer tablet, no waiting required.

To learn more about Indi with Snap + Core First, please visit tobiidynavox.com/indi



Indi shown with optional Durable Case.

Technical Support and Warranty

Service and support is easily available from one place, Tobii Dynavox, for all your Indi with Snap needs. Indi comes with a 2-year manufacturer's warranty and 90 days of telephone support to get you up and running quickly.



The optional Support360SM warranty comes with telephone support and covers any device repairs, no questions asked! Support360 is available for \$299 for one year, \$549 for two years and \$799 for three years. Alternatively, you are able to purchase telephone only support for \$99 per year.

Product Specifications

Type/Model	Tobii Dynavox Indi
Operating system	Microsoft® Windows 10 Pro
Screen	10.1" (25.65 cm)
Screen Resolution	1920 x 1200 px
Touch Panel	Multi-point Capacitive Touch
Weight	1 lb 11 oz (765 g)
Dimensions	9.6 x 6.9 x 1.4 in (24.4 x 17.6 x 3.5 cm)
Speakers	2 x 31mm 4 ohm, 3.5 W
Microphone	1 x Analog Microphone
Processor	Intel Atom Z8350 – 4 Cores, 2MB Cache, 1.92Ghz
RAM	4 GB RAM
Flash Storage	64 GB eMMC
Connectors	2 Micro-B USB 2.0 1 Micro-B USB 3.0 Headphone/Microphone Jack 12 VDC Charging Port

Optional Durable Case Accessories

Custom-built for AAC, an optional Durable Case is available in three different colors – Teal Blue, Raspberry, and Cool Gray – and comes in two different options, with or without keyguard holder.

- Durable Case - shoulder strap, and built-in mount plate, carry handle, and kickstand: \$99
- Durable Case with keyguard holder - shoulder strap, and built-in mount plate, carry handle, and kickstand: \$119
- Keyguards (for Durable Case with keyguard holder) are available for all different sizes of buttons: \$75

Infrared (IR)	Wideband Learning Infrared Remote Center
Sensors	3-Axis Gyro, Accelerometer, Compass, Ambient Light Sensor
Buttons	1 x Power On 1 x Volume Up and Volume Down
WLAN	IEEE 802.11 a/b/g/n
Bluetooth®	Bluetooth 4.0 / BLE
Camera	Rear-facing: 5 Megapixels Front-facing: 2 Megapixels
Power & Battery	Built-in 29.6 Wh Lithium Ion Battery
Battery Run Time	> 8 hrs. normal usage
Battery Charge Time	< 4 hrs
Windows Language Support	Supports all Windows 10 languages and language packs (download may be required)
Core First Language Support	English

Powerful built-in speakers

