

VISULAS YAG III Combi. Intelligent versatility.



The safe, compact way to combine retinal, cataract and glaucoma therapy.



Design + Flexibility = Office Convenience

VISULAS YAG III Combi is a big performer with a small footprint. Both laser systems can be placed safely and neatly under the instrument table, providing the user with more space and room for maneuver. The system's versatility is an added benefit. One example: While the VISULAS 532s is being used in the OR, the VISULAS YAG III can continue to be used independently. Both laser systems can be used as separate systems at any time.

During therapy, the vertically adjustable IT 4L laser table provides outstanding comfort for all patients, regardless of their size. Easy wheelchair access is guaranteed.

Simplicity + Intuition = Ease of Use

Both laser systems are reliably and safely controlled from a central console. Intuitive, language-independent user guidance facilitates the physician's work. Warnings and system information appear in clear text. The white on blue display is easily legible in a darkened room, allowing users with prescriptions of +/- 2.25 D to work without their glasses.

1 + 1 = 3

Add a photocoagulation laser to a disruption laser and what do you get? VISULAS YAG III Combi. The new, compact laser slit lamp from Carl Zeiss Meditec combines the benefits of the two systems and offers a spectrum of applications covering three areas of therapy.


VISULAS YAG III Combi is the optimum system not only for photocoagulation and post-cataract treatments, but also for such established glaucoma therapies as laser trabeculoplasty¹ and combined iridotomy². The push of a button is all you need to switch between the two laser systems.

VISULAS YAG III Combi is a three-in-one configuration: three therapies, one investment.



¹ Beckman H. et al; *The Glaucoma Laser Trial (GLT)-2, Results of Argon Laser Trabeculoplasty versus Topical Medicines*; *Ophthalmology* Nov 1990; 97, 1403-1413

² Lim L, Seah SK, Lim AS; *Comparison of argon laser iridotomy and sequential argon laser and Nd:YAG laser iridotomy in dark irises*; *Ophthalmic Surg Lasers* Apr 1996; 27(4), 285-288



Competence + Innovation = Technical Excellence

VISULAS YAG III Combi was specially developed to provide maximum safety and minimum fatigue. The optimized system for two laser sources features top-quality slit lamp optics from Carl Zeiss Meditec with the following superb benefits:

- The correct, highly transparent physician's safety filter is always active for all applications.
- The coaxial illumination and the servo-electric micromanipulator ensure that the coagulation laser is always visible in the center of the illuminated area to be treated.
- The parfocal system always images the laser spot sharply on the target area, ensuring homogenous, reproducible laser burns.
- The 4-point aiming beam and variable focus shift guarantee maximum precision in YAG laser therapy.
- Symmetrical arrangement of the instrument controls ensures ergonomic access to all slit lamp functions.

Needless to say, VISULAS YAG III Combi can also be used with our extensive range of proven accessories, e.g.

- Laser indirect ophthalmoscope
- Endo-probes
- DigCam adapter
- Coobservation tube and documentation devices





When your patients entrust you with their eyesight, their vision and your expertise converge. For the world's most advanced surgical and diagnostic solutions in ophthalmology, you can turn to Carl Zeiss Meditec. We're committed to earning your trust anew, every day.

Publication No.: 000000-1315-843
 Subject to change in design and scope of delivery and as a result of ongoing technical development.
 Printed in Germany, w.a.l 07/04. Printed on environmentally-friendly, unbleached paper.



Carl Zeiss Meditec AG
 Goeschwitzer Str. 51-52
 07745 Jena
 Germany

Phone: +49 (0) 36 41 / 22 0-3 33
 Fax: +49 (0) 36 41 / 22 0-2 82
 info@meditec.zeiss.com
 www.meditec.zeiss.com

Carl Zeiss Meditec Inc.
 5160 Hacienda Drive
 Dublin, CA 94568
 USA

Toll Free: (800) 342-9821
 Phone: +1 (925) 557-4651
 Fax: +1 (925) 557-4217
 info@meditec.zeiss.com
 www.meditec.zeiss.com