

Breast + Thorax Positioning

for accuracy, reproducibility + patient comfort



CIVCO[®]
Radiotherapy
a better way

MT-350

The MT-350 Breastboard combines reproducible breast positioning with efficient, comfortable patient setups. Its lightweight construction, including carbon fiber, provides ease-of-use with repeatable support.

MT-350 Features

- Grid treatment panels provide rigid, low-attenuating support
- Positive locking system offering angles in 2.5° increments from 0 to 20° and an angle of 25°
- Indexable hip stop prevents patients from sliding
- Wide variety of arm supports and headrest accessories
- Compatible with Standard, Extended and Tapered Wing Boards™



MT-350	MT350N
Wall Storage Device	MT350WSD



MT-350 Accessories - Sold separately

A	Standard Wing Board	MTWB09	E	Gel Head Cup	MT250SUB10
B	Tapered Wing Board w/ U-Grip (shown) w/ T-Grip	MTWB10U MTWB10T	F	Round Head Cup	MT250SUB11
C	Extended Wing Board w/ U-Grip w/ T-Grip (shown)	MTWB08U MTWB08T	G	Delrin Hand Grips	MT2504
D	Adapter Plate for Timo Headrests	MT350400	H	Bi-Axial Arm Support System™ Right Arm Support Left Arm Support Wrist Support	MTTAS2R MTTAS2L MTTWS2S