

Corrective Pitch Capability & Positioning Flexibility

Solstice™
SRS Immobilization System

FDA 510(k) Cleared!

Solstice™ SRS Immobilization System

- Variable pitch capability provides increased setup options in CT simulation and corrective positioning flexibility during treatment setup
- Midline thermoplastic mask attachment provides increased rigidity and minimizes shrinkage to enhance patient comfort
- Highly conformal to the patient when using a customizable cushion and dedicated thermoplastic mask
- Unique thermoplastic frame provides easy to use, secure and quiet mask attachment
- Utilizes ClearVision™ Open-face thermoplastic compatible with Vision RT's Align RT® camera system and IMRT Reinforced Thermoplastics™
- Adaptable with Precise Bite™ to aid in effective and repeatable positioning



Pitch at 0 (chin down)



Pitch at 5 (neutral)



Pitch at 10 (chin up)



Pitch indexing options 0 to 10

- Vision RT is unable to comment on, nor make claims, about the immobilization rigidity or patient alignment reproducibility when using this mask;
- Compatibility is subject to users following CIVCO guidelines for forming open masks in conjunction with Surface Guided Radiotherapy Systems, and subject to users following Vision RT's guidelines for drawing Regions of Interest for use in SRS with Open Masks;
- The mask vendor is responsible for ensuring that users are in receipt of guidelines for forming masks for use in conjunction with Surface Guided Radiotherapy Systems. Vendors may change their guidelines from time to time without reviewing these changes with Vision RT. Vision RT is not liable for misuse of masks;
- Vision RT makes no comment on, and is not liable for, suitability or otherwise of masks for SRS or any other procedure;
- Purchasers are responsible for making their own choices as to the masks selected by them.

CIVCO
Radiotherapy

System Customization

The system is comprised of a carbon fiber head support, customizable cushion and dedicated thermoplastic mask designed for secure and simple attachment. These components create a highly customizable and comfortable head and neck immobilization system for precise treatment to improve patient outcomes.

SGRT Starter Kit

visionrt COMPATIBLE

Kit includes: Solstice Head Support, Solstice ClearVision™ Open-face Thermoplastic, Standard (x10), 20 x 35 cm Chamfered AccuForm™ Cushion (x10), Precise Bite™ (x5).

CHS03SGK	Posifix®
CHS04SGK	Type-S™

Standard Starter Kit

Kit includes: Solstice Head Support, Solstice IMRT Reinforced Thermoplastic™ Style 18, Standard (x10), 20 x 35 cm Chamfered AccuForm Cushion (x10), Precise Bite (x5).

CHS03K	Posifix
CHS04K	Type-S

Individual Order Items

Head Support

CHS03	Carbon Fiber, Posifix (left image)
CHS04	Carbon Fiber, Type-S (right image)



ClearVision™ Thermoplastic

visionrt COMPATIBLE

MTAPCID1832SG	Standard (left), Qty 5
MTAPCIN1832SG	Narrow (right), Qty 5




IMRT Reinforced Thermoplastic™, Style 18

MTAPCID1832	Standard (left), Qty 5
MTAPCIN1832	Narrow (right), Qty 5



AccuForm™ Cushions

MTACL2035HF	20 x 35 cm Chamfered, Half Fill, Qty 1
MTACL2035	20 x 35 cm Chamfered, Qty 1



Precise Bite™

MTCBPBITE	Standard, Qty 5
MTCBPBITES	Small, Qty 5



Contact your sales representative for more information or to schedule a demo!



Global Sales Office

Orange City, Iowa, USA
800.842.8688 | +1 712.737.8688

info@CivcoRT.com
www.CivcoRT.com

COPYRIGHT © 2019. CIVCO, AND POSIFIX ARE REGISTERED TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. CLEARVISION, IMRT REINFORCED, PRECISE BITE, SOLSTICE AND TYPE-S ARE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. 2018P1357 REV. E

CIVCO[®]
Radiotherapy