

GE Healthcare

Lullaby* Resus Prime

Effective suction.
Controlled resuscitation.
One solution.

When it comes to infant resuscitation, time is as precious as your patients. Immediate access to all the right tools and equipment is essential.

The Lullaby Resus Prime supports:

- Oral/nasal pharyngeal suctioning as well as tracheal suction to patients.
- Monitored delivery of positive pressure through a manual resuscitation bag.†
- Clinical accuracy in pressure delivery.

This system is mountable on the Lullaby Warmer using a Dovetail-rail system and combines vital capabilities in a single, compact unit that is designed for high risk deliveries, intrahospital transport and bedside caregiving.

*Lullaby is a trademark of General Electric company
† Accessories like Airway manometer and PEEP valve may be required additionally

