

Manual operating table

For general and specialised surgery





PRACTICO

Wide variety of set-ups for different surgical procedures

The new hydraulic Practico is a versatile manually operated table model that enables a wide range of configurations for general and special surgery.

The operating table accommodates also larger patients, and offers excellent access to the operated site thanks to its design and adjustments options.

Additionally, the rust proof table top design makes the table suitable also for demanding climate conditions.

Practico manual can also be equipped with factory assembled Kidney Bridge seat/back section for kidney and thorax procedures. Reverse Trendelenburg



Features









Image: Contract of the second seco

Lots of benefits for the surgical team

- Excellent access for the surgical team to the patient
- Versatile use the same table for a wide variety of surgeries
- · Faster turnaround times thanks to easy set-up
- Compatible with many of the same accessories and table top sections as other Merivaara tables (Practico electro-mechanical and Promerix)
- C-arm and X-ray compatible
- Suitable for humid weather conditions thanks to corrosion-free material

Gas-spring assisted lateral tilt

Features

Divided leg sections (optional accessory)



Light foot rest for use during set-up (optional accessory)

Convenient kidney/thorax procedures



Hydraulic height adjustment

Patient positioning



Lateral position for kidney/thorax procedures with integrated kidney bridge seat/back section or with separate kidney elevator -accessory.



Lateral position for general surgery without kidney elevator.



Reverse mode for general or ENT surgery



Supine position for general surgery

PRACTICO

Patient positioning



Lithotomy position with stirrups and leg rests for, e.g., GYN and urology surgery



Lithotomy position with divided leg sections for, e.g., laparoscopic and bariatric surgery



Orthopedic extension device 19300 for trauma surgery

For further details on configurations and accessories, go to **www.merivaara.com** or contact your nearest Merivaara representative, see the website for contact info.



Easily removable sections thanks to the quick latching mechanism.



Available with VEF or IS mattresses

See **www.merivaara.com** for further details and possible updates.

- EN 60601-1, EN 60601-2-46 (when applicable)
- Practico manual is classified as class I medical device according to European Council; Directive 93/42/EEC Class I (Annex IX)

CE

Solutions for efficient physical patient flow and well-being

Merivaara provides a wide range of hospital-grade furniture such as opera-ting tables, medical lights, examination tables, trolleys and stretchers for transportation and day surgery, as well as birthing and patient beds. Additionally, our new concepts include both revolutionary integrated management systems and aftersales services.

Our versatile solutions and product offering enable convenient and cost-effective physical patient flow within the hospital, in day surgery and health clinics, as well as in nursing homes and home care. Our priority is to enhance comfort and well-being of both patients and their caregivers.

Established in 1901, Merivaara has more than 100 years of experience in designing and manufacturing hospital-grade furniture. Today

our products are highly appreciated by users in more than 120 countries thanks to their ease-of-use, durability and ergonomic design.

Merivaara medical products are CE-labelled and the company's quality management system complies with EU directives for medical device. It is certified by DNV according to ISO 9001:2008 and ISO 13485:2003 standards. The company's environmental management system is ISO 14001:2004 certified.



Merivaara Corp. Puustellintie 2, FI-15150 Lahti, FINLAND Tel. +358 3 3394 611, fax +358 3 3394 6144 merivaara@merivaara.com www.merivaara.com

The status of the information, specifications and illustrations in this publication is indicated by the date of publication. Whilst the greatest care has been taken over its compilation, Merivaara Corp. is not responsible for any errors or omissions. Merivaara Corp. reserves the right to make changes to the technology, features, specifications and design of the equipment and models without notice. Please check the Merivaara website for the latest updated specifications and other information.

All trademarks are Merivaara property unless otherwise stated. Patents granted or pending.

©Merivaara Corp. 6008EN 012013, printed in markprint 41580 2000



