

LOCK IN. ELEVATE YOUR CARE.

I-DEC™ Duodenoscope ED34-i10cT2

Reduces the risk of cross-contamination
without compromising maneuverability
or image quality.



Setting the new standard for ERCP with the only Disposable Elevator Cap (DEC™)



A **reliable therapeutic performance** is supported by a combination of an improved design of the **elevator shape** with the **ergonomic control** of a reusable endoscope.


Expect **outstanding visualisation** supported by a brighter HD+ image and finely differentiated colours and contrasts for proper positioning **to facilitate easier cannulation**.

A **reduced risk of cross-contamination** is ensured through a **disposable elevator cap**, DEC™ and increased surface access for cleaning of the distal end. This can reduce the distal end reprocessing steps up to 35% for more efficiency.



Inspired by #PENTAXMedicalTripleAim

SPECIFICATIONS	I-DEC™ DUODENOSCOPE
Model	ED34-i10cT2
Description	PENTAX Medical Video Duodenoscope ED34-i10cT2
Direction of view	102° (Retroflexed View 12°)
Field of view	100°
Depth of field	4–60 mm
Tip angulation up – down	120°/90°
Tip deflection right – left	105°/90°
Distal end width	13.6 mm
Insertion tube width	11.6 mm
Minimum instrument channel width	4.2 mm
Insertion tube working length	1250 mm
Compatible distal end cap	OE-A63 DEC Disposable Elevator Cap
Compatible processors	EPK-i8020c INSPIRA, EPK-i5500c IMAGINA



**EMEA Headquarter
Germany**

PENTAX Europe GmbH
Julius-Vosseler-Straße 104, 22527 Hamburg
T +49 (0)40 561 92-0 • F +49 (0)40 561 92-247
info.emea@pentaxmedical.com • www.pentaxmedical.com



**Manufacturer
HOYA Corporation**

6-10-1 Nishi-shinjuku
Shinjuku-ku, Tokyo
160-0023 Japan

PENTAX®
MEDICAL

CE Medical device class IIa – This product must be used only by healthcare professionals. Before use and for detailed product specifications, please refer to the instructions. In the interest of technical process, specifications may change without notice. PENTAX Medical INSPIRA™ and DEC™ are trademarks of HOYA Corporation. © Copyright 2023 HOYA Corporation. All rights reserved.

MKGI-3890-EN-G Rev. 3 (v3.1)