

Furniture suite

Sleeper Chair, Bedside Stand and Overbed Table

Ergonomic Sleeper Chair

MAB-510



Net weight	Width	Chair length	Bed length	Chair height	Bed height
52kg	74cm	100cm	200cm	96cm	58cm

Key features

- Easily rolls from armchair to bed position
- Two-step bed-to-chair folding system
- Integrated storage bin beneath surface
- Fully covered with vinyl
- Available in 24 different colors
- Easy-to-clean and durable design

Sleeper Chair

MAB-505



Net weight	Width	Chair length	Bed length	Chair height	Bed height
47kg	78cm	87cm	193cm	85cm	60cm

Key features

- Fully covered with vinyl
- Available in 24 different colors
- Easy-to-clean and durable design

Sleeper Chairs are available in 24 colors:



Your complete furniture solution

Two different types of Sleeper Chairs, both fully covered with easy-to-clean vinyl. Available in 24 different colors.

Bedside Stands include cross-lockable casters.

Overbed Tables contain dual locking casters and a gas spring to adjust the height of the table.

In combination with the SV1 bed, the furniture suite offers a complete and contemporary hospital room.

Bedside Stand

MAC-308 S



Key features

MAC-308 S

- Modern style with two drawers and mid shelf
- Steel body supported with aluminum profile
- HPL material for top surface and front panel
- Both drawers with removable ABS plastic shelving tray for enhanced cleanability
- Shelving tray dividers for ease of organization
- Cross-lockable casters for stability
- Push handle for ease of moving
- Available in five colors

MAC-307 S



Key features

MAC-307 S

- Modern style with single drawer and cabinet
- Steel body supported with aluminum profile
- HPL material for top surface and front panel
- Drawer with removable ABS plastic shelving tray for enhanced cleanability
- Shelving tray dividers for ease of organization
- Cross-lockable casters for stability
- Push handle for ease of moving
- Available in five colors

MAC-305-3



Key features

MAC-305 3

- Traditional style with 3 drawers
- Laminated chipboard body supported with aluminum profile
- ABS plastic top surface with raised edge for spill containment
- Cross-lockable casters for stability
- Available in five colors

Net weight	Width	Length	Height
29kg	46cm	47cm	81cm

Net weight	Width	Length	Height
25kg	46cm	47cm	81cm

Net weight	Width	Length	Height
31kg	52cm	46cm	83.5cm

Overbed Table

MOT-406



Key features

MOT-406

- Single-Slanting, Swivel Top with C-Shaped Base
- Compact HPL top surface
- Gas spring mechanism for height adjustability
- Dual locking casters for stability
- Available in five colors, with choice of gray or white frame

Net weight	Top width	Base width	Top length	Total length	Height
17kg	40cm	47cm	77cm	81cm	80.5-105cm

MOT-404



Key features

MOT-404

- Fixed Top with U-Shaped Base
- ABS plastic surface with raised edge for spill containment
- Gas spring mechanism for height adjustability
- Dual locking casters for stability
- Available in standard gray color only

Net weight	Top width	Base width	Top length	Total length	Height
15kg	42cm	34cm	80cm	83cm	80-104cm

All Bedside Stands and MOT-406 Overbed Table are available in 5 colors:

